

To aid with the assignment of signals in ^{13}C NMR spectra “Polarization Transfer” experiments can be performed to find out how many hydrogen atoms a carbon atom is carrying. These experiments begin with exciting all protons, followed by a transfer of the magnetization to the carbon atoms. Hence, only carbon atoms that bear protons will appear in these spectra.

1. The INEPT experiment

INEPT (Insensitive Nuclei Enhanced by Polarization Transfer) is no longer used because of frequently occurring phase distortions and has been replaced by the DEPT experiment.

2. The DEPT Experiment


DEPT (Distortionless Enhancement by Polarization Transfer) delivers clean ^{13}C subspectra which are easily phased, and which allow a reliable identification of CH, CH_2 , and CH_3 carbons. Quaternary Carbon atoms will not be seen in DEPT spectra.

Three experiments can be performed, which are identified by the last proton pulse angle:

DEPT-135 : CH: positive CH₂: negative CH₃: positive

DEPT-90 : CH: positive

~~DEPT-45 : CH: positive CH₂: positive CH₃: positive~~ ← **unnecessary**




Sample: 50 mg cholesteryl acetate in 0.5 ml C_6D_6

A combination of DEPT-135 and DEPT-90 will unambiguously reveal the CH multiplicities.

3. The Attached Proton Test (APT)

The APT experiment uses Spin-Echo instead of Polarization Transfer to produce spectra like DEPT-135. In addition, quaternary carbon atoms will appear as negative peaks (see arrows).



4. Which experiment(s) should I choose?

The best experiment(s) for your sample will depend on your exact situation. Some scenarios:

1. CH₂ present? DEPT-135 will identify CH₂ carbons as negative peaks
DEPT-90 will show only CH signals,
the remaining DEPT-135 signals will be CH₃ carbons.
2. No CH₂? APT will identify the quaternary carbons,
DEPT-90 will show only CH signals,
the remaining APT signals will be CH₃ carbons.


DEPT-45 is unnecessary. Instead, run DEPT-135 with more scans for a better sensitivity.

I highly recommend a two-dimensional experiment (HSQC). See the next page.

5. Heteronuclear Single Quantum Correlation (HSQC)

The HSQC experiment is much more sensitive than DEPT and is often the only option to get ^{13}C information from samples that are too dilute for direct ^{13}C observation.

HSQC spectra are color-coded and will allow you to distinguish between CH_2 groups (red crosspeaks) and CH and CH_3 groups (blue crosspeaks). Additionally, HSQC crosspeaks will show unambiguously which proton is attached to which carbon.



The arrows indicate the correlations for one CH (blue) and one CH_2 group (red).