


NMR Tips - Topspin Processing shortcuts

Bruker's Topspin has several commands that may help you speed up routine processing tasks. Some of these commands will only work (properly) in Topspin version 4.1 and higher.

Automatic phase and baseline correction:

apbk

^1H , ^{13}C : improvement of poor phase correction
 ^{11}B , ^{19}F , ^{29}Si : flattening of broad background signals and baseline distortions




Integration alternative:

presently for ^1H NMR only

sigreg

Before starting a manual integration, try:



Test sample: 10 mM cellobiose octaacetate in CDCl_3

SIGREG will only work with proton spectra.

Peak Picking:**pp**

Type “pp” and change one or both of the following two parameters to determine the number of peaks to be picked:

Minimum Intensity MI:

- larger values will skip smaller peaks.
- smaller values will find smaller peaks.
(in the noise between signals)

Detection Sensitivity PC:

- larger values will skip smaller peaks.
- smaller values will find smaller peaks.
(within each multiplet)

Required parameters	
Left picking limit F1P =	7.3188
Right picking limit F2P =	7.1226
Intensity of reference peak CY [rel] =	15
Minimum intensity MI [rel] =	0.1
Maximum intensity MAXI [rel] =	10000
Detection sensitivity PC =	10
Fraction of peak height for width calc. [0...1] =	0.5
Pick peaks of sign PSIGN =	pos. ▾
Reference peak selection mode PSCAL =	sreg ▾
Region file for PSCAL = sreg/psreg: SREGLST =	1H.CDCI3 ▾

First optimize **MI** to avoid picking empty regions. (= relative height of the smallest interesting peak)
Then optimize **PC** to adjust the number of peaks that are being picked above the MI level.

Calibration of any spectrum other than ¹H:**xiref**

Electronic chemical shift referencing is done mathematically with the IUPAC Xi scale. This command requires an accurately calibrated proton spectrum of the same sample. It also requires that all the spectra must be recorded in the same session as the ¹H spectrum.

1. Calibrate your ¹H spectrum using the signal of the solvent or of TMS (ideally).
2. Open your uncalibrated spectrum. Any experiment and any nucleus will work, i.e., COSY, NOESY, HSQC, HMBC, ¹³C, ¹⁹F, or any other nucleus.
3. Type the command “**xiref 10**” (assuming your ¹H spectrum is #10).

XIREF will allow you to accurately calibrate ¹³C NMR spectra if your solvent is D₂O. No additional internal or external standard is needed.